


This guide will walk you through configuring your iOS or iPadOS™ device with the VETSCAN VUE Rapid Test Analyzer. Before setup, you must install the VETSCAN VUE app and activate your VETSCAN VUE account. **If you have not received a welcome email with a link to activate your account, please contact Zoetis Technical Support at 1 800 822 2947 or dxsupport@zoetis.com before continuing.**

1 Log in to your VETSCAN VUE account

Open the VETSCAN VUE app. Press **Skip** to access the login screen. Enter the email you registered with Zoetis and the password you set up when you activated your VUE account. Press **LOGIN**.

Please note: Your mobile device must be connected to the internet to log in to the VETSCAN VUE app and complete setup.



2 Turn on your VETSCAN VUE

Follow the on-screen prompts to start device setup. Select **NEW DEVICE**.

Unbox your VETSCAN VUE, connect it to power, and turn it on.

Once the color light on your VETSCAN VUE stops flashing, press **CONTINUE**.



3 Connect to your VETSCAN VUE Rapid Test Analyzer

Follow the prompts to connect to the VETSCAN VUE on your mobile device.


Open **SETTINGS** on your mobile device and select **Wi-Fi**.

Select the network ID VUExxxx that matches AP: VUExxxx on the underside of your VETSCAN VUE.

Enter the AP password shown on the bottom of your VETSCAN VUE, then press **JOIN**.

Confirm that you are connected to the Wi-Fi network, VUExxxx, then return to the VETSCAN VUE app.


Press **CONTINUE** in the app.



4 Connect your VUE to the internet

Choose how to connect your VETSCAN VUE to the internet in your clinic. Press **CONNECT TO CORP. WI-FI** to use a Wi-Fi connection.

Or press **CONNECT TO ETHERNET** to use an ethernet connection.



Use a Wi-Fi connection

Enter the Network ID and password of your clinic's Wi-Fi. Press **CONNECT**.

Please note: Both the network name and password are case-sensitive. You may verify the password you enter using the Show Password slider.


Connecting may take several seconds. While connecting, the indicator light on the VUE will change colors. When the indicator light is solid green, press **CONTINUE**.

Use an Ethernet Connection

Follow the prompts to connect your VUE to your network using an ethernet cable, then press **CONTINUE** and **DONE**.

Connecting may take several seconds. While connecting, the indicator light on the VUE will change colors. When the indicator light is solid purple, press **CONTINUE**.

-OR-



5 Connect your mobile device to Wi-Fi


Follow the prompts to connect your mobile device to your clinic's Wi-Fi.

Open **SETTINGS** on your mobile device and select **Wi-Fi**.

Select your clinic's network ID. If prompted, enter the network password and press **JOIN**.

Confirm that you are connected to your clinic's Wi-Fi network, then return to the VETSCAN VUE app.

Press **CONTINUE** in the app.



6 Confirm that your VETSCAN VUE is connected

Confirm that the color light on your VETSCAN VUE is a solid green (Wi-Fi connection) or purple (Ethernet connection).

Press **FINISH SETUP**.



Congratulations, you are now ready to scan a VETSCAN Rapid Test with your VUE!